INFRARED FOREHEAD THERMOMETER user's manual

TABLE OF CONTENTS

Product Features2
 Product use2
 Technical Performance3
 Product description4
 How to use 6
 Principle of infrared temperature measurement10
 Body temperature common sense 11
Product Warranty / Warranty Card.....15

* Thank you for purchasing infrared forehead thermometer.
 * In order to use this product safely, please use it before use .Be sure to read the instruction manual.
 * After reading, please keep it in a safe place so that you can refer to it at any time.
 * Comes with a warranty card, do not lose it.
 * This product is powered by 3V, please use alkaline batteries.

Model: DY-A200
Jiangxi Dayou
Medical Technology Co., Ltd

1/24

Instructions

Precautions:
 -Follow the maintenance recommendations in this manual.
 -The ambient temperature of the product is + 10 C ~+40 C.
 -The product must be kept clean and placed in a dry place square.
 -Do not place the electronic thermometer in a place prone to electric shock square.
 -Do not place the electronic thermometer in an extreme temperature ring Environment: Above 50 C or below -20 C.
 -Do not place the electronic thermometer in a place where the humidity is higher than 80% RH surroundings.
 -Do not touch the lens inside the tip with your fingers.
 -Do not expose the product to sunlight or immerse it in water.
 -Do not drop the product.
 -If you find any problems, please contact the seller,Repair products.
 -Do not use in an environment with electromagnetic interference.
 -Please handle the product life according to local laws and regulations end-of-life waste and residue.

6/24

Features:

Accuracy: The measurement deviation is <0.3 C.
 Fast: Measurement time <1 second.
 Easy to use: One-button measurement, easy to operate.
 Non-contact: Measure on forehead without touching human skin and avoid cross infection.
 Measuring distance: 3cm ~5cm
 Large screen display: large screen LCD display.

High temperature alarm: freely set alarm temperature. Setting modification: Modify the setting parameters to suit different skins Colored race (white, black, yellowPeople, etc.).
 Unit conversion: Celsius, Celsius, Fahrenheit change.

Product Usage:
 Human body temperature measurement: measure human body temperature.
 Skin temperature measurement: Measure the surface temperature of human skin. Object temperature measurement: measure the surface temperature of an object, such as The temperature of the tea cup.
 Liquid temperature measurement: measure the temperature of liquids, such as baby wash
 Temperature of bath water, liquid in bottle
 Temperature etc.

2/24

Battery installation:

Use 2 batteries, pay attention when installing The positive and negative terminals of the battery cannot be reversed.

Use correctly:

The correct method of use is to test the accuracy the key is, otherwise it may cause measurement error.
 Because the infrared temperature measurement requires the surrounding environment high, so please follow the tips below.
 1. When measuring body temperature, point the instrument at the forehead Center — above the center of the eyebrow and stay vertical, measure the measuring area must not be covered by hair, the instrument and the forehead the recommended distance is about 3cm ~ 5cm.

2. When the person to be measured comes from a place with a large difference in temperature from the measurement environment, they should stay in the test environment for at least 5 minutes, and measure after the temperature is consistent with the

7/24

Technical performance

- Normal use conditions:
 Ambient temperature: + 10C~+40C Relative humidity: <85% RH.
- Power supply:
 DC 3V
- Product size:
 149x77x43mm(length X width X height)
- Product weight: 107g (excluding battery)
- Measurement range:
 Body temperature mode 32 C ~42.9 C Surface mode: 0 C ~ 100C
- Accuracy: 32 C ~ 34.9 C = ±0.3 C
 35 C ~ 42 C = ±0.2 C
 42.1 C ~42.9 C = ±0.3 C

- Power consumption: <50mW
- Measuring distance: 3cm ~5cm
- Automatic shutdown: is seconds. Temperature offset set to
 Press "set" button for 3 seconds, the screen showed that F1, press "set" button twice to enter the F3, then press "+" key to increase the temperature value of 0.1 °C (0. 1 °F), press key to reduce the temperature value of 0.1 (0. 1 °F). If the season or the environment changes, the forehead thermometer should be checked and adjusted.

Instructions

- environment, otherwise the measurement results will be affected.
- Fever patients' cold foreheads, sweating, and pluckingtaking other cooling measures will make the measurement results low should avoid measurement in this case.
 - The difference between the temperature of the instrument and the environment to be measured when taking out a larger place, it should be used placing the instrument in the environment of use 20Reuse in minutes.
 - The environment around the person being tested should be stable, can' t be at the air outlet of fan or air conditioner wait for the test where the air flow is heavy.
 - Not in a sunny place use the instrument.
 - It is recommended to measure 3 times each time when measuring right to show the most the data shall prevail.

8/24

Product description

4/24

Instructions

Alarm temperature adjustment: This product can adjust the default alarm temperature to suit the needs of different users. In the factory has done the default alarm temperature factory settings, factory default alarm temperature of37.5degrees C, if not necessary, it is recommended not to modify the factory default alarm value, if there is a need to modify, please follow the steps to adjust.
 Temperature alarm point setting Press "set" button for 3 seconds, the screen shows F1, press "set" button once to enter F2, then press key to increase the temperature value of 0.1 P (0. 1 °F), press key to reduce the temperature value of 0.1 °C (0. 1 °F), note: the default value of alarm is 37.5P .
 Temperature mode settings: Switching of body temperature or surface temperature patterns are divided into body temperature mode and surface temperature mode, which is used to measure body temperature, and surface temperature mode is used to measure the surface temperature of an object.
 Press the "set up" button to switch the temperature measurement mode.

9/24

Product description

5/24

Instructions

Temperature mode measurement range: 32 degrees C to 42.9 degrees C Surface temperature mode measurement range: 0 Cto100C Important: The surface temperature of the body watch is the temperature within the body is not the same, to measure the internal temperature must make use body temperature mode. Surface temperature pattern values are suitable for measurement external temperature, (e.g. bottles, bathrooms, indoor temperature Degrees, etc...)

Battery replacement:
 Products using 2 section battery, theoretically can be used more than 20,000 times in a row, when the screen has a battery symbol, and flash, indicating that the battery is low, need to replace the battery as soon as possible. Note:
 1. Open the battery cover to replace the battery, especially be careful that the battery is the negative pole should be placed correctly, otherwise there may not be Function.

10/24

Instructions

Precautions:
 -Follow the maintenance recommendations in this manual.
 -This product is suitable for professional use or home use.
 -The ambient temperature of the product is + 10 C ~+40 C.
 -The product must be kept clean and placed in a dry place square.
 -Do not place the electronic thermometer in a place prone to electric shock square.
 -Do not place the electronic thermometer in an extreme temperature ring Environment: Above 50 C or below -20 C.
 -Do not place the electronic thermometer in a place where the humidity is higher than 80% RH surroundings.
 -Do not touch the lens inside the tip with your fingers.
 -Do not expose the product to sunlight or immerse it in water.
 -Do not drop the product.
 -If you find any problems, please contact the seller,Repair products.
 -Do not use in an environment with electromagnetic interference.
 -Please handle the product life according to local laws and regulations end-of-life waste and residue.

6/24

Instructions

- When the product does not work, the power consumption is very small, when not used for a long time, it is recommended remove the battery to prevent the battery leak from damaging the instrument.Daily maintenance during the use of the product, follow the prompts when you find the following.
 - External dirt: wipe the dirt with clean soft cloth sticky water or wipe with cotton swab swabs, with a doctor wipe with alcohol can also also have the effect of sterilization, pay attention not to water or alcohol too much to avoid flowing into the meter damage inside the device.
 - Internal dirt: the black glass lenses of the internal probe are heavy to device, do not touch or top pressure, otherwise it will affect the accuracy of the measurement, when you find glass, when the lens surface is dirty, stick it with a cotton swab.95% of waterless alcohol wipes the surface of the lens.

11/24

Instructions

Important: You cannot wipe the lenses with 75% of the sterilizing alcohol (there will be traces of residual water). Lenses cannot be wiped with other chemical agents (which can cause damage to lenses).

Keeping and storing

Place in a dry, dark, sunny place.

Infrared temperature measurement principle

Understanding the principle of infrared temperature measurement can help you correctly with this instrument, the test data is more accurate.

- All objects, solids, liquids or gases are to the surrounding ring radiation releases infrared energy.
- The temperature of the object is proportional to the intensity of the radiated energy, i.e. the temperature the higher the amount of radiation, the greater the radiation energy.
- The energy of the human body's outward radiation is basically based on infrared radiationMain.

12/24

Body temperature common sense Measurement method Normal temperature

Anal temperature	Oral	36.6°C ~38°C
temperature	Underarm	35.5C ~37.5°C
temperature	Cochlear	34.7°C ~37.3°C
temperature	Forehead	35.8°C ~38°C

Changes in body temperature

Humans belong to the constant temperature animals, body temperature is basically constant, but not static, the body temperature is constantly changing in a day, as follows:

Evening Minimum due to sleep, reduced activity, lowest body temperature (low at 37 degrees C);
Morning, gradually high, from the warm bedding to the low room temperature in the morning, full the body muscles contract to produce heat; Noon high ester after meal, body temperature is the highest, the body naturally regulates; Night minimum due to the sun down, room temperature drop;

19/24

Body temperature common sense

- The infrared intensity of the human body's radiation to the outside world and the body itself body temperature is related to the ambient temperature of the human body, the higher the body temperature,tfte greater the radiation energy, the lower the ambient temperature, the less radiant energy bigger.
- By measuring the strength of the infrared energy radiated from the human body to the surrounding area the temperature of the human body itself can be calculated.
- The infrared thermometer itself does not emit any energy outwardsjust passively perceive the infrared energy of the outside world.
- infrared thermometer by accurate measurement of human release the faint infrared energy that comes out, and then through the complex transport calculation and various compensation corrections, can accurately measure the person body temperature. Products are made of a built-in infrared probe and phase the hardware, software composition, can be perceived, analyzed and remembered recorded objects and

13/24

Body temperature common sense Female

basic body temperature Women's body temperature and men are different, generally higher than men about 0.3 degrees C, ovulation period body temperature will be higher than usual 0.3 degrees C to 0.5 degrees C, so use this to know whether women are in ovulation, in order to facilitate fertility or contraception.

Definition: Basal Body Temperature, also known as resting body temperature, is the temperature measured by a woman after six to hours of sleep, such as waking up from a sleep in the morning, unaffected by exercise diet or mood changes. The base body temperature is usually the lowest body temperature in the body during the day and night.

The correct method for basic body temperature measurement

- Set up a thermometer or infrared thermometer to master correctly how to use it, for accuracy.
- Place the thermometer or infrared thermometer on the thermometer or THE infrared thermometer every night before going to bed wake up within reach of the place.

20/24

Body temperature common sense ambient

temperature. Therefore, once the user will close to a specific part of the human body (forehead) and press the test measure key, you can immediately activate the infrared radiation sensor, speed through the passive infrared sensor to detect due to arterial blood heat generated by the flow, so that the human body's Temperature.

Infrared temperature meter operation is simple, measurement time is only 1 second, easy to use, safe, not easy to bring the human body cross-infection, with infrared temperature meter, measuring temperature is simpler.

The main reasons for the screen to show Lo:

Why Lo information is displayed	Suggestions
Hair or sweat when the temperature is read	Make sure there are no obstacles on the forehead
There's cold air blowing the forehead.	Make sure no cold air blows your forehead directly
The forehead just colded.	Just cold applied for 10 minutes before measuring
Hair or sweat when the temperature is read	Make sure there are no obstacles on the forehead

14/24

Body temperature common sense

- After waking up in the morning every morning, put the thermometer under your tongue for 5 minutes to take out to read, the IR thermometer is aimed at the head of the forehead, and recorded on a specially crafted form.
- Before measuring body temperature, it is strictly forbidden to get up, stool, eat, talk, etc.
- It should be recorded whether there are many factors that affect the basic body temperature, such as: cold, Insomnia, drinking alcohol, taking medication, etc.

How to determine pregnancy by basic body temperature

The base body temperature usually rises with ovulation, the first half of the menstrual cycle, the body temperature fluctuates below 36.6 degrees C, and after ovulation is transferred to the menstrual cycle in the second half of the year, the body temperature is higher than the previous half, and fluctuates from 36.6 to 37 degrees C.Between, the increase in body temperature indicates ovulation, which lasts about two weeks after a period of high temperature, the

21/24

Body temperature common sense

The human body is a very complex biological integrated system. Body temperature is one of the most important data for whether human life activities are normal. Here are some medical temperature measurement common sense to help us understand our health status at any time.

Body temperature

Definition: The internal temperature refers to the true temperature at the pulmonary artery or esophagus of an organ in the body. Measurement method: Use a dedicated device with a temperature probe to measure deep into the pulmonary artery or esophagus. Advantages: Precise, constant, the most realistic response to the temperature of human internal organs.

Disadvantages: This is an immersed temperature measurement method that requires special equipment and expertise.

Anal temperature

Definition: Refers to the temperature of the human rectum.

15/24

Body temperature common sense body

temperature drops again, and then it is the turn of the monthly menstruation. If pregnant, the base body temperature does not drop, but continues to be in high temperatures, which are due to an increase in the amount of hormones secreted by the ovaries after pregnancy. But a fever in a woman with poor ovaries or a cold can affect the underlying body temperature curve, so the accuracy is not 100 percent.

Attached:The property description, ownership change and interpretation rights of this product belong to Jiangxi Da-you Medical Technology Co., Ltd

22/24

Body temperature common sense

Measurement method:

Use a thermometer to measure the anorectum. Advantages: The anal temperature is relatively stable and close to the internal temperature of the human body.generally speaking, the body temperature refers to the anal temperature. Disadvantages: inserted into the rectum, not easy to accept, easy to cross infection.

Oral temperature

Definition: Refers to the temperature under the tongue of the human mouth. Measurement method: using a thermometer Insert under the tongue of the mouth and close your mouth for measurement. Advantages: The oral temperature is relatively stable, close to the internal temperature of the human body. Disadvantages: inserted into the mouth, not easy to be accepted, easy to cross-infection.

Underarm temperature

Definition: Refers to the temperature at the armpit spans of the human body.

16/24

17/24

PRODUCT GUARANTEE

- The product will be free for one year with the purchase invoice from the date of purchase maintenance.
 - We will notFree maintenance is available.
- (A) Unauthorized disassembly and modification of products caused by the failure;
(B) Failure caused by accidental fall during use and handling;
(C) Failure caused by lack of reasonable maintenance;
(D) Failure to follow the correct instructions in the instructions for usefault;
(E) Failure caused by improper repairs in repair shops not authorized by our company wait a minute.
- Repair services outside the scope of maintenance will be charged in accordance with the regulations.
- When requesting maintenance services, please consult the after-sales service department.

23/24

Measurement method:

Measured under the armpit with a thermometer. Pros: Simple measurement, not easy to cross infection, easy to be accepted, is the most common temperature measurement methods are used in clinical practice. Cons: Lower than actual body temperature, each measurement changes measurements are not accurate enough.

Cochlear temperature

Definition: Refers to the temperature of the inner eardrum of a human cochlea. Measurement method: The cochlear implant measurement is inserted with an infrared ear temperature thermometer. Advantages: The temperature of the eardrum in the cochlea is relatively stable, close to the human body temperature. Disadvantages: Infrared temperature measurement method is required, susceptible to the size of the ear canal and the ear the effects of the taorad (commonly known as ear poop), resulting in measurement deviations, after measuring 1 to 2 times, you must wait 3 to 5 minutes before you can

17/24

18/24

CERTIFICATE	
Product Model: _____	Name: _____
Date of purchase: ____	Body code: _____
Address: _____	
_____ phone: _____	
Sales store name: _____	

Registrant/Manufacturer:

Jiangxi Dayou Medical Technology Co., Ltd

Address: Building B7,Yangguang Avenue
Standard Workshop,Wanzai industrial Park,
Wanzai County, Yichun City, Jiangxi
Province.

Production License No.:
Jiangxi Food and Drug Administration Machinery
Production Permit No. 20200315
Registration Certificate No.:
Jiangxi Machinery Injection 20202070074 Medical
Device Product Technical Requirement No.: Jiangxi
Machinery Injection 20202070074 After-sales service
unit:

Jiangxi Dayou Medical Technology Co., Ltd Address:
Building B7, Yangguang Avenue
Standard Workshop, Wanzai industrial Park,
Wanzai County, Yichun City, Jiangxi
Province.
Service hotline: 0752-2488887

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Reason for displaying "LO" or "HI"	suggestion
Impact by hair or sweat	Make sure no block when testing
Impact by the unstable air	Make sure the air is stable
Impact by the far distance	Pay attention to the distance(1cm-3cm)
From outdoors to indoors with big difference of temperature	Please wait for 20 minutes after the human body fit environment temperature

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amfori BSCI
Trade with purpose

Product Features

- Non-Contact and with a high accuracy of measurement
- Two kinds of measurement mode selection: body temperature and surface temperature
- Exchange between °C and °F
- Function of sound notification of high body temperature
- Capable of storing 32 sets for the latest measurement data
- Saving the data and auto turn off
- Function of auto selection the range
- LCD display resolution :0.1 °C(0.1°F)
- Backlight LCD digital display

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